



Anti-GFP

Product Name	Anti-GFP
Catalog No	IT-000-010
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	Green Fluorescent Protein (GFP)
Applications	Western Blotting (1:100-1:1000); Immuno-precipitation
Specificity	Reacts with Green Fluorescent Protein. Cross-reacts to most GFP derivatives
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

