



# Anti-Luciferase (*Metridia*)

---

Product Name	Anti-Secreted luciferase ( <i>Metridia</i> )
Catalog No	IT-000-012
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	<i>Metridia longa</i> secreted luciferase, a 24 kDa protein from the marine copepod <i>Metridia longa</i> (Genebank No. AAR17541)
Applications	Western Blotting (1:200-1:1000)
Specificity	Reacts <i>Metridia longa</i> secreted luciferase protein. Cross-reactivity not tested.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

