

# Anti-Cypridina Luciferase



Product Name	Anti-Cypridina Luciferase
Catalog No	IT-000-013
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	Recombinant Cypridina Luciferase (a.a. 1-168) protein (Genebank No. BAD08210)
Applications	WB(1:100), ELISA(1:1000-1:100,000), IP, etc.
Specificity	Reacts with Cypridina Luciferase protein. Minimal cross-reactivity to other luciferase proteins.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

