

Anti-Cypridina Luciferase



Product Name	Anti-Cypridina Luciferase
Catalog No	IT-000-013M4
Source	Mouse monoclonal antibody
Isotype	IgG1
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol
Immunogen	Recombinant Cypridina Luciferase (a.a. 1-168) (Genebank No. BAD08210)
Applications	WB, ELISA, IP, IF etc.
Epitope	Linear
Specificity	Reacts with Cypridina Luciferase. No Cross-reactivity to other luciferase protein.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing.</u> The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

ELISA Profile of IT-000-013M4

