



Anti-Gaussia Luciferase

Product Name	Anti-Gaussia Luciferase
Catalog No	IT-000-014
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	Recombinant Gaussia Luciferase (a.a. 1-185) (Genebank No. AAG54095)
Applications	WB(1:1000), ELISA(1:1000-1:100,000), IP, etc.
Specificity	Reacts with Gaussia Luciferase protein. Minimal cross-reactivity to other luciferase proteins.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing</u> . The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

