

# Anti-*Gaussia* Luciferase



Product Name	Anti- <i>Gaussia</i> Luciferase
Catalog No	IT-000-014M1
Source	Mouse monoclonal antibody
Isotype	IgM
Formulation	Each vial contains 100 µg purified IgM in 100 µl PBS with 40% glycerol
Immunogen	Recombinant <i>Gaussia</i> Luciferase (a.a. 1-168) (Genebank No. AAG54095)
Applications	ELISA, IP, IF etc.
Epitope	Conformational
Specificity	Reacts with <i>Gaussia</i> Luciferase. No Cross-reactivity to other luciferase protein.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

## ELISA Profile of IT-000-014M1

