



Anti-6xHis Tag

Product Name	Anti- 6xHis Tag
Catalog No	IT-000-His-M2
Source	Mouse monoclonal antibody, clone 121E11
Isotype	IgG1
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol
Immunogen	KHL-6xHis
Applications	WB, ELISA, IP, IF etc.
Epitope	Linear
Specificity	Reacts with 6xHis Tag in various tagged recombinant proteins. Suitable for Western Blotting, ELISA, IP and IF, etc.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

