

Anti-Human IgG



| | |
|----------------|--|
| Product Name | Anti-Human IgG |
| Catalog No | IT-000-hIg-M3-HRP |
| Source | Mouse monoclonal antibody, clone 2H6 |
| Isotype | IgG2a |
| Formulation | Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol. |
| Immunogen | DNA vaccine expressing human IgG1 Fc protein |
| Applications | ELISA, IP, IF etc. |
| Epitope | Conformational |
| Specificity | HRP labeled monoclonal antibody. Reacts with human IgG and monkey IgG. Minimal cross-reactivity to IgG from others species as tested. |
| Storage | Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required. |
| Use Limitation | For research use only, not for use in diagnostic procedures. |