



Anti-gp120/160(HIV-1/Clade A)

Product Name	Anti-gp120/160 (HIV-1/Clade A)
Catalog No	IT-001-001
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	gp140 Env protein (a.a. 1-486) (HIV-1/Clade A/) (Genebank No.AF084657)
Applications	Western Blotting (1:200-1:1000)
Specificity	Reacts with HIV-1(Clade A) gp120/160 proteins. Cross-reactivity to envelope protein from other clades not tested.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

