



# Anti-gp120/gp160(Pan)

Product Name	Anti- gp120/gp160(Pan)
Catalog No	IT-001-002M27
Source	Mouse monoclonal antibody, clone 6B3
Isotype	IgG2a
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	Recombinant gp120 (SF162)(B) protein (IT-001-0028p) (a.a. 30-503) (Genebank No. <a href="#">AAT67508</a> ).
Applications	WB, ELISA, IP, IF etc.
Epitope	Linear
Specificity	Reacts with most HIV-1 gp120/gp160 proteins as tested. Minimal cross-reactivity to other subtypes.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

## ELISA Profile of IT-001-002M27

