




# gp120 (HIV-1/Clade C)

Product Name	gp120 (HIV-1/Clade C)
Catalog No	IT-001-003p
Source	Viral protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS only or PBS containing 0.1% BSA and 25% glycerol.
Protein	6xHis tagged HIV-1 gp120 protein (a.a. 24-515) of strain C.ZA.1197MB from South Africa (Genebank #AY463234)
Applications	Western Blotting standard, antibody ELISA, immunogen, etc.
Purity	≥ 95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: purified gp120 (Clade C) protein.