



# gp41(Ectodomain)(HIV-1)

---

Product Name	gp41( $\Delta$ TM)(Bal)
Catalog No	IT-001-0051p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 $\mu$ g of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged gp41 (Bal)(HIV-1/Clade B) protein (a.a. 537-676) (Genebank No. <a href="#">M68893</a> )
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing</u> . The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

