




# Anti-gp41(HIV-1/Clade C)

Product Name	Anti-gp41(Clade C)
Catalog No	IT-001-006
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	HIV-1 gp41(a.a.512-684) (Genebank No. AY463234)
Applications	Western Blotting (1:200-1:1000)
Specificity	Reacts with gp41(Clade C) protein. Cross-reactivity to other clades not tested.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



C: 293 cell extract control  
E: 293 cell expressing gp41(Clade C)

Note: