

Anti-p24 (HIV-1)



| | |
|----------------|--|
| Product Name | Anti-p24(HIV-1) |
| Catalog No | IT-001-007M9 |
| Source | Mouse monoclonal antibody, clone 40G10 |
| Isotype | IgG2a |
| Formulation | Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol. |
| Immunogen | DNA vaccine expressing Gag protein (Clade A) of HIV-1(Genebank No. AY772985) |
| Applications | WB, ELISA, IP, IF etc. |
| Epitope | Linear |
| Specificity | Reacts with p24 and Gag proteins from HIV-1. |
| Storage | Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required. |
| Use Limitation | For research use only, not for use in diagnostic procedures. |

