




Anti-Gag(Pan)(HIV-1)

Product Name	Anti-Gag(Pan)
Catalog No	IT-001-007PAN
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	HIV-1 Gag(Con_Con) (http://www.hiv.lanl.gov/cgi-bin/CONSENSUS_DOWNLOAD/ConsensusDownloader.cgi).
Applications	Western Blotting (1:200-1:1000)
Specificity	Reacts with HIV-1 Gag protein. Recognizes Gag polypeptide from Clades A, B, C HIV-1 viruses.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



C: 293 cell extract control
1: 293 cell expressing Gag (Clade C)
2: 293 cell expressing Gag (Clade B)
3: 293 cell expressing Gag (Clade A)

Note: