

# Anti-p24 (HIV-1)



Product Name	Anti-p24(HIV-1)
Catalog No	IT-001-009M8
Source	Mouse monoclonal antibody, clone 57B7
Isotype	IgG1
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	DNA vaccine expressing Gag protein (Clade C) of HIV-1(Genebank No. AAB61122)
Applications	WB, ELISA, IP, IF etc.
Epitope	Linear
Specificity	Reacts with p24 and Gag proteins from HIV-1.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

## ELISA Profile of IT-001-009M8

