




p17 (HIV-1/Clade A)

Product Name	p17 (HIV-1/Clade A)
Catalog No	IT-001-010Ep
Source	Viral protein purified from E.coli
Formulation	Each vial contains 100 µg of lyophilized protein (1 mg/ml) in PBS with 8M urea. Reconstitute the protein with 100µl of Millipore water
Protein	C-terminal 6xHis tagged p17 (HIV-1/Clade A) matrix protein (Genebank accession no.AF004885)
Applications	WB standard, antibody ELISA, etc
Purity	>90%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: p17 (HIV-1/Clade A) protein

Note: