




# Anti-Nef (HIV-1/Clade A)

Product Name	Anti-Nef (HIV-1/Clade A)
Catalog No	IT-001-013
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	Nef (HIV-1/Clade A) protein (Genebank accession #AF084657)
Applications	Western Blotting (1:200-1:1000)
Specificity	Reacts with Nef protein from HIV-1/Clade A viruses. Cross-reactivity to Nef from other clades not tested.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



C: 293 cell extract control;  
E: 293 cell expressing Nef

Note: