




Anti-Nef (HIV-1/CladeB)

Product Name	Anti- Nef (Clade B)
Catalog No	IT-001-014
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	HIV-1 Nef(B) (HXBc2) (Genebank No. K03455)
Applications	Western Blotting (1:200-1:1000)
Specificity	Reacts with HIV-1 Nef(B). Cross-reactivity to Nef from other clades not tested.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



C: 293 cell extract control
E: 293 cell expressing Nef(B)

Note: