



Anti-Nef(Pan)(HIV-1)

Product Name	Anti-Nef (Pan)(HIV-1)
Catalog No	IT-001-014PAN
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	Nef protein(Con-Con) (http://www.hiv.lanl.gov/cgi-bin/CONSENSUS_DOWNLOAD/ConsensusDownloader.cgi)
Applications	Western Blotting (1:200-1:1000), ELISA(1:1000-5000)
Specificity	Reacts with HIV-1 Nef protein from clade A, B, C.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

