

Anti- gp120(SIV)(mac239)



Product Name	Anti- gp120(SIV)(mac239)
Catalog No	IT-001-022-M1-HRP
Source	Mouse monoclonal antibody, clone 36-E2
Isotype	IgG1
Formulation	Each vial contains 50 µg HRP labeled IgG in 100 µl PBS with 40% glycerol
Immunogen	Recombinant SIV gp120 protein (mac239)
Applications	ELISA, WB, IP, IF etc.
Epitope	Linear
Specificity	Reacts with SIV gp120(mac239)

Storage Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.

Use Limitation For research use only, not for use in diagnostic procedures.

