



Gag(SIV-1/mac239)

Product Name	Gag(SIV-1/mac239)
Catalog No	IT-001-024Ep
Source	Viral protein purified from E.coli
Formulation	Each vial contains 100 µg of protein (1 mg/ml) in PBS with 8M urea.
Protein	C-terminal 6xHis tagged SIV-1 Gag(mac239) protein (Genebank No. M33262)
Applications	WB standard, antibody ELISA, etc
Purity	>90%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing</u> . The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

