

Anti-gp120/gp160(A244)(Clade 01_AE)



Product Name	Anti-gp120/gp160(A244)(Clade 01_AE)
Catalog No	IT-001-165M2
Source	Mouse monoclonal antibody, clone 26A4
Isotype	IgG1
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	Recombinant gp140(A244)(Clade 01_AE) protein (a.a. 44-681)(IT-001-165p)
Applications	WB, ELISA, IP, IF etc.
Epitope	Linear
Specificity	Reacts with gp120/gp140/gp160(A244)(Clade 01_AE), and some Clade A gp120 proteins as tested. Minimal cross-reactivity to other subtypes.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

ELISA Profile of IT-001-165M2

