

Anti-gp120 (HIV-2/B_UC-1)



Product Name	Anti- gp120 (HIV-2/b_UC-1)
Catalog No	IT-001-203M4
Source	Mouse monoclonal antibody, clone 7E3
Isotype	IgM
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	DNA vaccine expressing gp120 (HIV-2/B_UC1) protein (a.a.22-513) (Genebank #: Q76638).
Applications	WB, ELISA, IP, IF etc.
Epitope	Linear
Specificity	Reacts with gp120 (HIV-2/B_UC-1) only. Not reactive to gp120 (HIV-2/NIH-2) and gp120 (HIV-2/Ben). Cross-reactivity to others not tested.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

ELISA Profile of IT-001-203M4

