



# Biotinylated hCD4(Ectodomain)

---

Product Name	Biotinylated hCD4(Ectodomain)
Catalog No	IT-001-CD4p-Biotin
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 50 µg of biotinylated protein (0.5 mg/ml) in PBS+40% glycerol
Protein	6xHis tagged human CD4 ectodomain (a.a.1-398) (GenBank No. P01730).
Applications	antibody ELISA, protein labeling, etc
Purity	>95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.