




# Anti-Spike (SARS-CoV)

Product Name	Anti-S (SARS-CoV)
Catalog No	IT-002-001
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	Spike protein of SARS-CoV (Beijing02) (Genebank No. AF372624)
Applications	Western Blotting (1:100-1:1000)
Specificity	Reacts with spike protein of SARS-CoV. Cross-reacts to most spike proteins from other subtypes of SARS-CoV.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



C: 293 cell extract control  
E: 293 cell expressing S(Beijing02)

Note: