

Anti-S1(MERS-CoV/EMC 2C/2012)



Product Name	Anti-S1(MERS-CoV/EMC 2C/2012)
Catalog No	IT-002-010M7
Source	Mouse monoclonal antibody, clone 8B2
Isotype	IgG1
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	DNA vaccine expressing S1 of the spike protein (MERS-CoV/EMC 2C/2012)(Genebank No. AFS88936)(a.a. 18-725)(IT-002-010p).
Applications	ELISA, IP, IF etc.
Epitope	Conformational
Specificity	Reacts with spike protein (MERS-CoV/EMC 2C/2012) only as tested.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

