



Anti-S1(Bat coronavirus HKU4-2)

| | |
|----------------|--|
| Product Name | Anti-S1(Bat coronavirus HKU4-2) |
| Catalog No | IT-002-020M3 |
| Source | Mouse monoclonal antibody, clone 20F9 |
| Isotype | IgG1 |
| Formulation | Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol. |
| Immunogen | DNA vaccine expressing S1 of the spike protein (Bat coronavirus HKU4-2)(Genebank No. ABN10848)(a.a. 21-729)(IT-002-020p). |
| Applications | ELISA, IP, IF etc. |
| Epitope | Linear |
| Specificity | Reacts with spike protein (Bat coronavirus HKU4-2) only as tested. |
| Storage | Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required. |
| Use Limitation | For research use only, not for use in diagnostic procedures. |

