



Anti-S1(Bat coronavirus HKU4-2)

Product Name	Anti-S1(Bat coronavirus HKU4-2)
Catalog No	IT-002-020M4
Source	Mouse monoclonal antibody, clone 26A9
Isotype	IgG1
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	DNA vaccine expressing S1 of the spike protein (Bat coronavirus HKU4-2)(Genebank No. ABN10848)(a.a. 21-729)(IT-002-020p).
Applications	ELISA, IP, IF etc.
Epitope	Conformational
Specificity	Reacts with spike protein (Bat coronavirus HKU4-2), and weakly cross-reacts with spike protein (MERS-CoV/EMC 2C/2012).
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

