



SARS Core (340-390 a.a.)

Product Name	SARS Core (340-390 a.a.)
Catalog No	IT-002-052Ep
Source	Recombinant protein purified from E.coli
Formulation	Each vial contains 100 µg of protein in PBS with 8M Urea.
Protein	The E.Coli derived 32kDa recombinant protein contains the Nucleocapsid core protein 340-390 amino acids immunodominant regions.
Applications	WB standard, antibody ELISA, immunogen, etc.
Purity	>95% , based on SDS PAGE
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

