




Anti-HA(H1N1)

Product Name	Anti-HA(H1N1)
Catalog No	IT-003-001
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	Hemagglutinin (amino acid 1-343)(H1N1)(A/New Caledonia/20/99) (Genebank No. ABF21272)
Applications	Western Blotting (1:200-1:1000)
Specificity	Reacts with human and avian H1(H1N1). Cross-reactivity to other HAs not tested.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing.</u> The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



C: 293 cell extract control
E: 293 cell expressing H1

Note: