


HA(Δ TM)(A/Christchurch/16/2010) (H1N1)



Product Name	HA(Δ TM)(A/Christchurch/16/2010)(H1N1)
Catalog No	IT-003-00103 Δ TMp
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 μ g of purified protein (1 mg/ml) in PBS
Protein	C-terminal 6xHis tagged HA(A/Christchurch/16/2010)(H1N1) (a.a.18-529) (See a.a. sequence below)
Applications	Antibody ELISA, Immunogen, WB standard, etc.
Purity	>95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: HA(Δ TM)(A/Christchurch/16/2010)(H1N1)

Protein Sequence:

DTLCIGYHANNSTDTVDTVLEKNVTVTHSVNLLEDKHNGKLCCLRGVAPLHLGKCNIAGWILGN
PECESLSTASSWSYIVETSSSDNGTCYPGNFIDYEELREQLSSVSSFERFEIFPKTSSWPDHDSNKG
TAACPHAGAKSFYKNLIWLVKKGNSYPTLSKSYINDKGKEVLVLWGIHPSTADQQSLYQNAD
AYVFGTSRYSKKFKPEIAIRPKVRNQEGRMNYYWTLVEPGDKITFEATGNLVAPRYAFAMERNA
GSGIISDTPVHDCNTTCQTPKGAINSLPFQNIHPITIGKCPKYVKSTKLRLATGLRNVPSIQSRGLF
GAIAGFIEGGWTGMVDGWYGYHHQNEQSGYAADLKSTQNAIDKITNKVNSVIEKMNTQFTAV
GKEFNHLEKRIENLNKKVDDGFLDIWTYNAELLVLENERLTDYHDSNVKNLYEKVRSQKLNNA
KEIGNGCFEFYHKCDNTCMESVKNNGTYDYPKYSEEAKLNREEIDGVKLESTRIYQ