




HA₁ (H1N1)(PR/8/34)

Product Name	HA ₁ (H1N1/PR8)
Catalog No	IT-003-0010p
Source	Viral protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS containing 0.1% BSA and 25% glycerol.
Protein	6xHis tagged hemagglutinin (amino acid 18-345) (H1N1) (A/PR/8/34/(H1N1)) (Genebank No. P03452)
Purity	≥ 95%
Applications	Western Blotting standard, antibody ELISA, immunogen, etc.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: purified HA₁ (H1N1/PR/8/34) protein.