


# HA( $\Delta$ TM)(H1N1)(A/Puerto Rico/8/1934)



---

Product Name	HA( $\Delta$ TM)(H1N1)(A/Puerto Rico/8/1934)
Catalog No	IT-003-0010 $\Delta$ TMp
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 $\mu$ g of purified protein (1 mg/ml) in PBS
Protein	C-terminal 6xHis tagged HA( $\Delta$ TM)(H1N1)(A/Puerto Rico/8/1934) protein (a.a.18-529) (Genbank # P03452)
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: HA( $\Delta$ TM)(H1N1)(A/Puerto Rico/8/1934)