


Anti-HA1(A/Beijing/262/95)(H1N1)



Product Name	Anti- HA1(A/Beijing/262/95)(H1N1)
Catalog No	IT-003-00110M2
Source	Mouse monoclonal antibody, clone 10E3
Isotype	IgG2b
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	DNA vaccine expressing HA1(A/Beijing/262/95)(H1N1) protein (a.a. 18-344)(Genebank accession AAP34323).
Applications	ELISA, IP, IF etc.
Epitope	Linear
Specificity	Reacts with H1 only. Minimal cross-reactivity to other subtypes.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

ELISA Profile of IT-003-00110M2



ELISA Profile of IT-003-00110M2

