

Anti-H1N1/H3N2(Seasonal)



| | |
|----------------|--|
| Product Name | Anti-H1N1/H3N2(Seasonal) |
| Catalog No | IT-003-001M11-HRP |
| Source | Mouse monoclonal antibody, clone 12B9 |
| Isotype | IgG1 |
| Formulation | Each vial contains 50 µg HRP labeled IgG in 100 µl PBS with 40% glycerol |
| Immunogen | Recombinant HA proteins from seasonal strains |
| Applications | WB, ELISA, IP, IF etc. |
| Epitope | Linear |
| Specificity | Reacts with most seasonal flu strains. No cross-reactivity to other subtypes. |
| Storage | Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required. |
| Use Limitation | For research use only, not for use in diagnostic procedures. |

