



Anti- HA(H1N1/H3N2(Seasonal))

Product Name	Anti-HAH1N1/H3N2(Seasonal)
Catalog No	IT-003-001M11
Source	Mouse monoclonal antibody, clone 12B9
Isotype	IgG1
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol
Immunogen	Sequential immunizations using recombinant H1(ΔTM) (H1N1/A/New Caledonia/20/99)(a.a.18-530)(Genebank No. AAP34324), and H5(H5N1)(A/Vietnam/1203/2004) protein (a.a. 18-530) (GeneBank No. AY651334)
Applications	WB, ELISA, IP, IF, etc.
Epitope	Linear
Specificity	Reacts with most seasonal flu strains. No cross-reactivity to other subtypes.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

ELISA Profile of IT-003-001M11

