



Anti- H1(Seasonal 1999-2010)

Product Name	Anti- H1(Seasonal 1999-2010)
Catalog No	IT-003-001M26
Source	Mouse monoclonal antibody, clone 18D7
Isotype	IgG2b
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	DNA vaccine expressing HA1 (A/New Caledonia/20/99(H1N1)) protein (a.a. 18-345) (IT-003-001p).
Applications	WB, ELISA, IP, IF etc.
Epitope	Linear
Specificity	Reacts with seasonal H1 (1999-2010) as tested. Cross-reacts with H9. Minimal cross-reactivity to other subtypes.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing.</u> The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

