


# Anti-H3(A/Victoria/361/2011)(H3N2)



Product Name	Anti- H3 (A/Victoria/361/2011)(H3N2)
Catalog No	IT-003-00423M12
Source	Mouse monoclonal antibody, clone 26B3
Isotype	IgG2b
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	DNA vaccine expressing HA1(A/Victoria/361/2011)(H3N2) protein (a.a. 18-348)(IT-003-00423p).
Applications	ELISA, IP, IF etc.
Epitope	Conformational
Specificity	Pan H3 as tested. Minimal cross-reactivity to other subtypes.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

**ELISA Profile of IT-003-00423M12**



**ELISA Profile of IT-003-00423M12**

