

Anti-H3(H3N2)(Seasonal)



Product Name	Anti-H3(H3N2)(Seasonal)
Catalog No	IT-003-004M2
Source	Mouse monoclonal antibody, clone 111E9
Isotype	IgG1
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol
Immunogen	Recombinant HA proteins (H3N2)(A/Brisbane/10/2007)
Applications	ELISA, WB, IP, IF etc.
Epitope	Linear
Specificity	Reacts with recent seasonal flu H3N2 (Wyo, Wis, Bris, Victoria) strains (also see the datasheet of IT-003-004M1). No cross-reactivity to other subtypes.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

