

## Anti-H3(H3N2)( (A/Brisbane/10/2007))

---

Product Name	Anti-H3(H3N2)(A/Brisbane/10/2007)
Catalog No	IT-003-004M3
Source	Mouse monoclonal antibody, clone 64D1
Isotype	IgG2b
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol
Immunogen	Recombinant HA proteins (H3N2)(A/Brisbane/10/2007)
Applications	ELISA, Neut., IP, IF etc.
Epitope	Conformational
Specificity	Reacts with H3(H3N2)(A/Brisbane/10/2007). Weakly reacts with H3 (H3N2)(A/Wisconsin/67/X-161/2005). No cross-reactivity to other subtypes.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

