



Anti-HA1(A_Ontario_RV1273_2005)(H3N2)

Product Name	Anti- HA1(A_Ontario_RV1273_2005)(H3N2)
Catalog No	IT-003-004SW5M2
Source	Mouse monoclonal antibody, clone 71D5
Isotype	IgG2b
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	Recombinant HA1(A_Ontario_RV1273_2005)(H3N2) protein(Gene Bank Accession No. ABF17954)(human isolate of SOIV H3N2) (IT-003-004SW2p)
Applications	WB, ELISA, IP, IF etc.
Epitope	Linear
Specificity	Reacts with swine origin human H3N2. Cross-reactive with some seasonal H3N2. Also cross-reactive with H4, H7, H14. Epitope in the head region.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

