

Anti-H5(A/Hubei/1/2010)(H5N1)



Product Name	Anti- H5(A/Hubei/1/2010)(H5N1)
Catalog No	IT-003-00519M10
Source	Mouse monoclonal antibody, clone 77E7
Isotype	IgG2b
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	DNA vaccine expressing H5 (A/Hubei/1/2010)(H5N1) protein (Genbank No. EPI267032) (IT-003-00519ΔTMp).
Applications	ELISA, IP, IF etc.
Epitope	Conformational
Specificity	Reacts with Pan H5, also some H7, H8 and H14 as tested. No Cross-reactivity to other subtype as tested.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

