

Anti- H5(Pan)(A/H5N1)



Product Name	Anti- H5(Pan)(A/H5N1)
Catalog No	IT-003-005M24
Source	Mouse monoclonal antibody, clone 18C7
Isotype	IgG1
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	DNA vaccine expressing HA1 (A/Vietnam/1194/2004(H5N1)) protein (a.a. 1-345) (IT-003-005p).
Applications	ELISA, IP, IF etc.
Epitope	Conformational
Specificity	Reacts with most H5 (H5N1) tested. Cross-reacts with some H7, H8 and H14. Minimal cross-reactivity to other subtypes.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

ELISA Profile of IT-003-005M24

