

Anti-H5



| | |
|----------------|--|
| Product Name | Anti- H5 |
| Catalog No | IT-003-005Pan-M4 |
| Source | Mouse monoclonal antibody, clone 2G5G6 |
| Isotype | IgG1 |
| Formulation | Each vial contains 100µg purified IgG in 100µl PBS with 40% glycerol |
| Immunogen | DNA vaccine expressing HA(ΔTM)(A/Indonesia/5/2005(H5N1)) protein (IT-003-0052ΔTMp). |
| Applications | WB, ELISA, IP, IF etc. |
| Epitope | Linear |
| Specificity | Reacts with most H5 as tested, not react with other subtype strains as tested. |
| Storage | Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required. |
| Use Limitation | For research use only, not for use in diagnostic procedures. |

