

Anti- HA(H9N2)



Product Name	Anti- HA(H9N2)
Catalog No	IT-003-009M2-HRP
Source	Mouse monoclonal antibody
Isotype	IgG1
Formulation	Each vial contains 50 µg HRP labeled IgG in 100 µl PBS with 40% glycerol
Immunogen	Recombinant HA1 protein (H9N2)(A/Guinea fowl/Hong Kong/WF10/99 (H9N2)) protein (a.a.20-339) (Genebank accession AAO46082).
Applications	ELISA, WB, IP, IF etc.
Epitope	Linear
Specificity	Reacts with most HA(H9N2) strains. No cross-reactivity to other subtypes.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

