



Anti-H13(H13N8)

Product Name	Anti-H13(H13N8)
Catalog No	IT-003-016M5
Source	Mouse monoclonal antibody, clone 7E1E5
Isotype	IgG2a
Formulation	Each vial contains 100µg purified IgG in 100µl PBS with 40% glycerol
Immunogen	Recombinant H13 protein (A/black-headed gull/Netherlands/1/00 (H13N8)) (IT-003-014ΔTMp)
Applications	ELISA, IP, IF etc.
Epitope	Conformational
Specificity	Reacts with H13 protein (A/black-headed gull/Netherlands/1/00 (H13N8)). The epitope is located within HA1 region. No cross-reactivity to other subtypes.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

