




# NP(H3N2)(A/Brisbane/10/2007)

Product Name	NP(H3N2)(A/Brisbane/10/2007)
Catalog No	IT-003-024Ep
Source	Viral protein purified from E.coli
Formulation	Each vial contains 100 µg of lyophilized protein (1 mg/ml) in PBS. Reconstitute the protein with 100µl of Millipore water.
Protein	C-terminal 6xHis tagged NP (H3N2)(A/Brisbane/10/2007) protein (Genebank accession No. <a href="#">EU199248</a> )
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: NP (H3N2)(A/Brisbane/10/2007) protein

Note: