


HA1(B/Phuket/3073/2013)



Product Name	HA1(B/Phuket/3073/2013)
Catalog No	IT-003-B11p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS
Protein	C-terminal 6xHis tagged HA1(B/Phuket/3073/2013) protein (a.a. 16-362) (see amino acid sequence below)
Applications	Antibody ELISA, Immunogen, WB standard, etc.
Purity	>95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing</u> . The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



Protein SEQ

driectgitsnsphvvtatqgevnvtgviplttptksyfanlkgtrtrgkclpdelnctldldvalgrpmcvgttptsakasilhe
vrpvtsgcfpimhdrtkirqlpnllrgyekirlstqnvidaeakapggpyrlgtsgscpnatskigffatmawavpkdnyknt
npltvvevpyictgedqitvwgfhsdntqmkslygdsnpqkftssangvtthyvsqigdfpdqtedgglpqsgrivvdy
mmqkpgktgtivyqrgvllpqkvwcasgrskvikgspligeadclheeygglnskspyytgkhakaigncpiwvktpl
klangtkyrppakllkerg