




HA(Δ TM)(B/Phuket/3073/2013)

Product Name	HA (Δ TM)(B/Phuket/3073/2013)
Catalog No	IT-003-B11 Δ TMp
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 μ g of purified protein (1 mg/ml) in PBS.
Protein	6xHis tagged HA(Δ TM)(B/Phuket/3073/2013)(See below for aa sequence)
Applications	WB standard, antibody ELISA, etc
Purity	\geq 95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing</u> . The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: HA(Δ TM)(B/Phuket/3073/2013)

Note:

Protein Seq:

driectgitssnsphvvktatqgevnvtgviplttptksyfanlkgtrtrgkclcpdclnctldlvalgrpmcvgtppsakasilhe
vrpvtsgcfpimhdkirqlpnllrgyekirlstqnvdaekapggpyrlgtsgscpnatskigffatmawavpkdnyknt
npltvvepyictegedqitvwgfhdsnktqmkslygdsnpqkftssangvtthyvsqigdfpdqtedgglpqsgrivdy
mmqkpgkgtgivyqrgvllpqkvwcasgrskvikgspligeadclhekygglnskpyytgehakaigncpiwvktpl
klangtkyrppakllkergffgaiagfleggwegmiagwhgytshgahgvavaadlkstqeainkitknlnslselevknlq
rlsgamdelhneileldekvdldradtissqielavllsnegiinsedehllalerklkmlgpsavdigngcfetkhkcnqtcl
driaagtdagefslptfdslnitaaslnddgldnh